

## Product specification affigene® DNA extraction

**Code No** 71000942

**Intended use** Sample preparation for affigene® trender/tracer real-time PCR assays

**Total kit size** 46 determinations (corresponds to a complete affigene® trender/tracer kit for 48 determinations)

**Prev. of false neg. results** Internal Control introduced during sample preparation

**Regulatory status** CE labelled according to IVD Directive 98/79/EC.

**Sample matrix** Serum  
Plasma (ACD, EDTA)  
Whole blood (EDTA)  
Vesicle swabs  
Throat swabs  
Cerebro Spinal Fluid (CSF)  
Broncho Alveolar Lavage (BAL)

**Sample volume** 50 µL required for whole blood protocol  
200 µL required for cell-free sample protocols